# Human IgG ELISA kit

#### IBT Bioservices cat# 0101-001, lot# 2304002

#### Instructions for use

#### **1.** Purpose:

For the quantitative measurement of human IgG.

#### **2.** Reagents supplied:

Reagent supplied	Lot Number	Protein	Amount	Storage		
		Conc.		Temperature		
DPBS 1X for diluting	N/A	N/A	13 mL	4°C		
Capture Antibody						
Capture Antibody	03.15.2023	1.0 mg/mL	30 µL	short-term 4°C		
				long-term -20°C		
Assay Diluent	N/A	N/A	80 mL	4°C		
for blocking						
Standard, Samples,						
Secondary Antibody						
& Detection Reagent						
Human IgG Standard	02.09.2019-В	1.0 mg/mL	15 μL	short-term 4°C		
				long-term -20°C		
Secondary Antibody	01.26.2021	0.474	30 µL	short-term 4°C		
		mg/mL		long-term -20°C		
Detection Reagent	N/A	0.500	15 μL	4°C		
		mg/mL				
TMB one-step	N/A	N/A	13 mL	4°C		
substrate						



### **3.** Reagents required but not included in the kit:

- DPBS powder (Corning cat# 55-031-PB) or equivalent, stored at 4°C, for preparing ELISA Wash Buffer. For more information, please refer to page 6.
- TWEEN-20 (Sigma-Aldrich cat# P1379-500 mL) or equivalent, stored at ambient temperature, for preparing ELISA Wash Buffer. For more information, please refer to page 6.
- Deionized water

#### 4. Materials required but not included in the kit:

- MaxiSorp flat-bottom, polystyrene, 96-well plates (Thermo Scientific cat# 439454)
- Polypropylene TiterTubes, maximum volume for each tube = 1 mL (Bio-Rad cat# 223-9391) or equivalent, used to prepare standard and sample dilutions
- Microplate sealing film
- Polypropylene 15 mL and 50 mL conical tubes
- Reagent reservoirs
- Absorbent papers

### 5. Equipment required but not included in the kit:

- Automatic plate washer (example: BioTek ELx450)
- Plate reader with capability of measuring absorbance at 650 nm (example: Molecular Devices plate reader)
- Software for graphing the Standard as a 4PL curve and for calculating the Unknown Samples from the standard curve (example: Softmax software)
- Single-channel and multi-channel pipettes



#### 6. Assay Procedure:

- 1. Prepare Capture antibody solution
  - Briefly spin the **Capture Antibody (1.0 mg/mL)** vial and gently mix by pipetting up and down
  - Dilute 1:500 in DPBS 1X to target 2 µg/mL
  - Example: For one full plate, add 24 μL Capture antibody to 12 mL of DPBS 1X
- Add 100 μL/well of Capture antibody solution to the MaxiSorp plate. Cover plate using plate sealing film. Incubate covered plate overnight at 4°C.
- 3. The following day, equilibrate plate and Assay Diluent to ambient temperature for at least 15 min.
- 4. Empty contents from the plate and wash 3 times (each time 300  $\mu$ L/well) with Wash Buffer using an automatic plate washer or multi-channel pipette.
- 5. Add 250  $\mu$ L/well of Assay Diluent to block non-specific binding. Incubate for at least 1 hour at ambient temperature.
- During blocking step, prepare dilutions of STANDARD and UNKNOWN SAMPLES in TiterTubes. Each dilution should contain sufficient volume (>200 μL) for duplicate wells.



- a. STANDARD
  - Briefly spin the Standard (1.0 mg/mL) vial and gently mix by pipetting up and down
  - First dilution = 1:10 to target 100 µg/mL

Add 10  $\mu L$  STANDARD at 1.0 mg/mL to 90  $\mu L$  of Assay Diluent. This is 10  $\mu g/mL$  INTERMEDIATE.

• Second dilution = 1:10 to target 10 μg/mL

Add 50  $\mu$ L INTERMEDIATE at 100  $\mu$ g/mL to 450  $\mu$ L of Assay Diluent. This is the 10  $\mu$ g/mL starting concentration tested on ELISA plate.

Use a new pipet tip to gently mix by pipetting up and down.

- Serial 1:3.162-fold dilutions
  - Transfer 139 μL from the previous dilution to 300 μL Assay Diluent.
  - Discard pipet tip
  - Use a new pipet tip to gently mix by pipetting up and down. Change pipet tips between each dilution.
  - Repeat for subsequent dilutions

#### b. UNKNOWN TEST SAMPLES

Note: Mouse and rat serum at dilutions 1:100, 1:1,000 and 1:10,000 did not interfere with the recovery of the Human IgG standard (please see page 9). Other matrices or lower dilutions have not been tested.

- Prepare sample dilutions using Assay Diluent at dilution factors determined by the end-user
- Remember to mix well by pipetting solution up and down a few times and change pipet tips in between dilutions
- The followings are suggestions for sample dilution factors:



Unknown	<b>Dilution Factor</b>	
Test		
Sample		
Dilution 1	1:100	1:100 dilution
		Add 5 μL Unknown Sample
		to 495 μL of Assay Diluent
Dilution 2	1:1,000	1:10 of 100 dilution
		Add 30 μL Dilution 1
		to 270 μL of Assay Diluent
Dilution 3	1:10,000	1:10 of 1,000 dilution
		Add 30 μL Dilution 2
		to 270 μL of Assay Diluent

- 7. Empty contents from the plate and wash 3 times (each time 300  $\mu$ L/well) with Wash Buffer using automatic plate washer or multichannel pipette.
- Use multi-channel pipettor to transfer 100 μL/well of STANDARD or UNKNOWN SAMPLE dilutions from TiterTubes to duplicate wells in MaxiSorp plate. Change pipet tips appropriately to avoid cross-contamination. Cover plate with plate sealing film. Incubate for 2 hours at ambient temperature.



- 9. At the end of the 2-hour incubation step, prepare Secondary Antibody solution
  - Briefly spin the Secondary Antibody (0.474 mg/mL) vial and gently mix by pipetting up and down
  - Dilute 1:474 in Assay Diluent to target 1 µg/mL.
  - Example: For one full plate, add 25.4 μL of Secondary Antibody to 12 mL Assay Diluent
- 10. Empty contents from the plate and wash 3 times (each time 300  $\mu$ L/well) with Wash Buffer using automatic plate washer or multichannel pipette.
- 11. Add 100  $\mu$ L/well of Secondary Antibody solution to plate. Cover plate with plate sealing film. Incubate for 1 hour at ambient temperature.
- 12. At the end of the 1-hour incubation step, prepare Detection Reagent solution
  - Briefly spin the **Detection Reagent** vial and gently mix by pipetting up and down
  - Dilute Detection Reagent 1:4,000
  - Example: For one full plate, add 3 µL of Detection Reagent to 12 mL Assay Diluent
- 13. Empty contents from the plate and wash 3 times (each time 300  $\mu$ L/well) with Wash Buffer using automatic plate washer or multichannel pipette.
- 14. Add 100  $\mu$ L/well of Detection Reagent solution to plate. Cover plate with plate sealing film. Incubate for 1 hour at ambient temperature, shielded from light.
- 15. During this time, equilibrate TMB substrate to ambient temperature, shielded from light.



- 16. Empty contents from the plate and wash 3 times (each time 300  $\mu$ L/well) with Wash Buffer using automatic plate washer or multichannel pipette.
- 17. Add 100  $\mu$ L/well of TMB substrate. Incubate plate at ambient temperature, shielded from light. Start timer for 30 min color development.
- 18. Immediately following the 30 min color development, place plate in the plate reader programmed to shake the plate for 5 sec prior to end-point read at 650 nm wavelength.
- 19. Prepare a standard curve from the data produced from the serial dilutions with concentrations on the x axis (log scale) vs. absorbance on the y axis (linear). Interpolate the concentration of the Unknown Samples from the standard curve.

#### Notes regarding plate washing:

- ELISA Wash Buffer (1X DPBS + 0.05% TWEEN-20):
  - Dissolve one bottle of DPBS powder in deionized water to prepare 10 liters of 1X DPBS
  - Add 5 mL TWEEN-20 to 10 L of 1X DPBS
  - o Gently mix
- Use BioTek plate washer model ELx405, "COSTAR\_FLAT" program (Number of cycles: 3; Volume wash buffer: 300 μL/well).
- Empty the MaxiSorp plate's content into biohazard container and blot on paper towels
- Wash plate using "COSTAR\_FLAT" program
- Tap plate on paper towels to remove any residual liquid.
- Immediately add solution to the wells. Do not let the wells dry for extended time.



### 7. Example Template and Standard Curve

	EXAMPLE OF PLATE TEMPLATE											
	1	2	3	4	5	6	7	8	9	10	11	12
A	STD	STD	STD	STD	STD	STD	STD	STD	STD	STD	STD	Blank
B	10,000	3,162	1,000	316	100	31.6	10.0	3.2	1.0	0.32	0.10	
C	Spl 1	Spl 1	Spl 1	Spl 4	Spl 4	Spl 4	Spl 7	Spl 7	Spl 7	Spl 10	Spl 10	Spl 10
D	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3
E	Spl 2	Spl 2	Spl 2	Spl 5	Spl 5	Spl 5	Spl 8	Spl 8	Spl 8	Spl 11	Spl 11	Spl 11
F	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3
G	Spl 3	Spl 3	Spl 3	Spl 6	Spl 6	Spl 6	Spl 9	Spl 9	Spl 9	Spl 12	Spl 12	Spl 12
H	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3	Dil.1	Dil.2	Dil.3

#### **EXAMPLE OF STANDARD CURVE**



#### DATA ANALYSIS

Softmax software is used to calculate Human IgG (ng/mL) detected in Unknown Samples, based on the 4PL standard curve, tested on every plate, using the following equation:

$$X = C \star \left( \frac{A - Y}{Y - D} \right)^{(1/B)}$$

X = ng/mL Y = Absorbance Value (OD 650 nm)

A = Lower asymptote B = Slope

C = Inflection point

D = Upper asymptote



Version v1, page 8

	HulgG standard diluted in					HulgG standard diluted in				
HulgG	Assay	Naïve Mouse Serum			HulgG	Assay	Naïve	e Mouse Serum		
ng/mL	Diluent	1:100	1:1,000	1:10,000	ng/mL	Diluent	1:100	1:1,000	1:10,000	
10000.00	3.329	3.336	3.326	3.315	10000.00	100%	100%	100%	100%	
3162.28	3.322	3.325	3.346	3.327	3162.28	100%	100%	101%	100%	
1000.00	3.416	3.359	3.372	3.250	1000.00	100%	98%	99%	95%	
316.23	3.394	3.325	3.346	3.212	316.23	100%	98%	99%	95%	
100.00	3.233	3.134	3.103	3.043	100.00	100%	97%	96%	94%	
31.62	2.579	2.474	2.457	2.266	31.62	100%	96%	95%	88%	
10.00	1.456	1.335	1.337	1.301	10.00	100%	92%	92%	89%	
3.16	0.582	0.585	0.572	0.543	3.16	100%	101%	98%	93%	
1.00	0.281	0.308	0.268	0.274	1.00	100%	110%	95%	98%	
0.32	0.153	0.159	0.161	0.159	0.32	100%	104%	105%	104%	
0.10	0.120	0.129	0.131	0.135	0.10	100%	108%	109%	113%	
	Mean	OD650 of	duplicate	wells		Comparing OD650 to the standard				
						dilu	ited in As	ay Diluen	t	
	Hu	IgG standa	rd diluted	in		HulgG standard diluted in				
HulgG	Assay	Naï	ve Rat Ser	um	HulgG	Assay	Naïve Rat Serum			
ng/mL	Diluent	1:100	1:1,000	1:10,000	ng/mL	Diluent	1:100	1:1,000	1:10,000	
10000.00	3.182	3.177	3.220	3.146	10000.00	100%	100%	101%	99%	
3162.28	3.179	3.150	3.178	3.160	3162.28	100%	99%	100%	99%	
1000.00	3.255	3.175	3.211	3.116	1000.00	100%	98%	99%	96%	
316.23	3.197	3.157	3.228	3.067	316.23	100%	99%	101%	96%	
100.00	3.000	2.963	2.938	2.870	100.00	100%	99%	98%	96%	
31.62	2.468	2.399	2.354	2.243	31.62	100%	97%	95%	91%	
10.00	1.431	1.424	1.332	1.319	10.00	100%	100%	93%	92%	
3.16	0.623	0.641	0.616	0.599	3.16	100%	103%	99%	96%	
1.00	0.279	0.321	0.257	0.251	1.00	100%	115%	92%	90%	
0.32	0.161	0.190	0.170	0.155	0.32	100%	118%	106%	96%	
0.10	0.122	0.163	0.139	0.139	0.10	100%	134%	114%	114%	
	Mean	OD650 of	duplicate	wells		Comparing OD650 to the standard				
						dilu	ited in As	ay Diluen	t	

## 8. Recovery of Human IgG spiked in rodent serum

